

2 WATTS

Repeat 29 Hz F2-4

20 msec  
40 mJ  
1160 mW Av  
58.0 % Duty

MODE

SUPER PULSE

Frenum Vascular

Custom Setting	Watts - SP/CW - Mode - Air
CS 1 FRENUM VASCULAR	02-CW-Repeat F2-4-High
CS 2 HEMOSTASIS	02-CW-Repeat F1-4-High
CS 3 MUCOSA INCISION	04-CW-Repeat F2-2-High
CS 4 FASCIA INCISION	04-SP-Repeat F2-2-High
CS 5 MUSCLE INCISION	02-CW-Repeat F2-4-High

light Scalpel

4 WATTS

Repeat 29 Hz F2-2

10 msec  
40 mJ  
1160 mW Av  
29.0 % Duty

MODE

SUPER PULSE

Fascia

2 WATTS

Repeat 20 Hz F1-4

20 msec  
40 mJ  
800 mW Av  
40.0 % Duty

MODE

SUPER PULSE

Hemostasis

4 WATTS

Repeat 29 Hz F2-2

10 msec  
40 mJ  
1160 mW Av  
29.0 % Duty

MODE

SUPER PULSE

Mucosa

2 WATTS

Repeat 29 Hz F2-4

20 msec  
40 mJ  
1160 mW Av  
58.0 % Duty

MODE

SUPER PULSE

Muscle